




Suportul Vital Pediatric

Suport Vital Avansat



IN TIMPUL RCP

- RCP de calitate: frecventa, compresiuni, decompresiuni
- Planificati actiunile inainte de intreruperea RCP
- Administrati oxigen
- Acces vascular (intravenos, intraosos)
- Adrenalina la 3-5 minute
- Protejarea cailor aeriene si capnografie
- Continuati compresiunile toracice dupa protejarea cailor aeriene
- Corectati cauzele reversibile

CAUZE REVERSIBILE

- Hipoxia
- Hipovolemia
- Hipo-/hiperpotasemia/dezechilibre metabolice
- Hipotermia
- Pneumotorax in tensiune
- Toxice
- Tamponada cardiaca
- Tromboembolism